

Long-Term Effects of a Dynamic Ankle Foot Orthosis on a Patient with Cerebral Palsy Following Ischaemic Perinatal Stroke – a Case Study

Daniel Sabbagh¹, Jörg Fior¹, Ralf Gentz¹

¹FIOR & GENTZ Gesellschaft für Entwicklung und Vertrieb von orthopädietechnischen Systemen mbH, Lüneburg, Germany

Background and Aim

Ischaemic perinatal stroke (IPS) is a common cause of cerebral palsy (CP) in new-born infants. Those children develop neurological deficits which lead to gait disorders. Ankle foot orthoses (AFOs) are an important factor in immediate and continual gait improvement for children with CP. For structuring the orthotic treatment, the pathological gait of CP patients can easily be assessed by using the Amsterdam Gait Classification (fig. 1) [1]. According to this classification, the characteristics of an AFO can be defined individually for each gait type [2]. A well-structured orthotic treatment can assist physiotherapeutic approaches and create new cerebral connections by using the right motor impulses [3].

Research question: Are there long-term changes in kinematic and spatio-temporal gait parameters by wearing dynamic AFOs?

Materials and Methods

TYPES OF GAIT ACCORDING TO THE AMSTERDAM GAIT CLASSIFICATION					
TYPES OF GAIT	Type 1	Type 2	Type 3	Type 4	Type 5
KNEE	normal	hyperextended	hyperextended	flexed	flexed
FOOT CONTACT	complete	complete	incomplete	incomplete	complete

Fig. 1: Amsterdam Gait Classification – pathological gait of CP patients is classified in 5 gait types which can be easily assessed by evaluating knee position and foot contact in mid stance. Our case study patient is gait type 5 (flexed knee and complete foot contact). This gait type is named 'crouch gait' as well.




Fig. 2: a) dynamic AFO, b) patient with fitted dynamic AFO

A 15-year-old female patient with CP following ischaemic perinatal stroke (hemiparesis of left body-side) classified as gait type 5 of the Amsterdam Gait Classification (fig. 1) is treated with a dynamic AFO (carbon-composite, adjustable range of motion, defined pivot point; fig. 2a). The AFO is adjusted to the patient's individual normal posture (fig. 2b) [4]. On the day of delivery (DD) and at the follow-up (FU) after 3 months, 2-dimensional gait analyses were performed – each with and without AFO. Thereby, measurements of hip, knee and ankle angle plus step/stride length, velocity, cadence, and time of stance and swing were taken.

Results




Fig. 4: joint kinematics of hip, knee and ankle on day of delivery without AFO (DD w/o AFO), with a dynamic AFO (DD w/ AFO), at follow-up without (FU w/o AFO) and with AFO (FU w/ AFO) compared to a reference group (gray band)




Fig. 5: The bar graphs show differences between conditions on DD w/o AFO, on DD w/ AFO at FU w/o AFO and at FU w/ AFO for spatio-temporal gait parameters. The table shows mean values and standard deviation.

Measurements on DD showed increased step length (+5%), longer single-supported (+8%) and shorter double-supported phase (-19%) with AFO. Ankle plantar flexion diminished in swing until loading response by 7° and knee flexion in stance by 5°. FU showed increased stride length (+5%), velocity (+8%) and cadence (+2%) plus reduced hip (-10°) and knee flexion (-3°) in stance with and without AFO, compared to DD.

Discussion and Conclusion

On DD the AFO improved kinematic and spatio-temporal gait parameters. Increasing ankle dorsiflexion in swing and initial contact, enabling passive plantar flexion at loading response and decreasing peak ankle dorsiflexion in mid stance were immediate effects of the dynamic AFO (tab. 1). During three months of wearing the AFO, the CP patient had come closer to a physiological gait which shows in a less flexed hip, knee and ankle during stance with and without AFO at FU – especially in mid stance. The dynamic properties of the AFO improved spatio-temporal gait parameters even without AFO. An adaptation process is stated by these long-term changes in gait but the presence of new cerebral connections after IPS must still be proved.

Tab. 1: Immediate effects of the dynamic AFO at DD and FU. Long-term effects are shown in the differences between results of DD w/o AFO and FU w/o AFO. Immediate effects are improvements of ankle dorsiflexion in swing and peak ankle dorsiflexion. Long-term effects can be seen in most kinematic and spatio-temporal gait parameters.

	Day of delivery (DD)		Follow-up (FU)	
	immediate effects		immediate effects	long-term effects
Kinematics				
Hip extension	slight increase	-	-	increase
Knee extension	slight increase	-	-	slight increase
Ankle dorsiflexion (ISw - LR)	increase	increase	increase	slight increase
Ankle dorsiflexion (MSt - PSw)	-	slight decrease	slight decrease	decrease
Peak ankle dorsiflexion	decrease	decrease	decrease	decrease
Spatio-temporal				
Step length [m]	slight increase	slight increase	slight increase	increase
Stride length [m]	slight increase	-	-	increase
Velocity [m/s]	slight increase	slight increase	slight increase	increase
Cadence [steps/min]	-	increase	increase	increase
Single Support [s]	slight increase	decrease	decrease	slight increase
Double Support [s]	slight decrease	-	-	slight increase

Keywords: orthotics, cerebral palsy, ischaemic perinatal stroke, gait analysis, dynamic AFO, long-term effects

Literature

[1] Becher J. Journal Prosthetics Orthotics 2002; 14(4): 143-149.
[2] Gentz R et al. Orthopädie Technik 2012; 63(8): 35-41.

[3] Horst R. Motorisches Strategietraining PNF 2005; Thieme.
[4] Sabbagh D et al. Gait Posture 2015; 42(S1): S80.