



Fig. 1

Application: for padding orthoses

Material properties: exchangeable and soft padding material. It is antimicrobial and can be washed at 60°C. The padding material has a skin-friendly towelling coating on one side and a velour coating on the reverse side. Use the hook dots or the hook tape (see below) for a safely fixing of the padding material in an orthosis. When cutting the material with pinking shears, there are no frayed edges. This means that you get a functional edge which does not need to be linked.

xDRY Towelling Padding Material (fig. 1)

Article Number	Scope of Delivery	Length x Width x Thickness [mm]	Colour	Unit
PL3697-02/1	1 x fig. 1	1000 x 1170 x 4	black	sheet
PL3697-02/2	1 x fig. 1	1000 x 1170 x 4	blue	sheet
PL3697-02/3	1 x fig. 1	1000 x 1170 x 4	beige	sheet
PL3697-02/4	1 x fig. 1	1000 x 1170 x 4	grey	sheet



Fig. 2

To fix the padding material, the self-adhesive hook tape or hook dots can be ordered separately.

Application: for fixing the towelling padding material in the orthosis

Material properties (hook tape): 1mm thin micro hook tape

Material properties (hook dots): 1mm thin micro hook dots with a diameter of 22mm

Hook Tape, Self-adhesive, Black (fig. 2)

Article Number	Scope of Delivery	Length x Width x Thickness [mm]	Unit
KV4050-L1000	1 x fig. 2	1000 x 50 x 1	piece



Fig. 3

Hook Dots, Self-adhesive, Black (fig. 3)

Article Number	Scope of Delivery	Unit
KV0020-22/06	6 x fig. 3	set
KV0020-22/40	40 x fig. 3	set