



## CROSSROADS

with Elastic Outsole Made of EVA Material (brown)  
with Elastic vibram® Flat Sole (black)

Shoe Measurements		Use
pitch*	15mm	h-Cast
heel height**	25mm	Orthosis Configurator

Size	36-48	
Material	High quality leather – certified cow leather used for the upper material and calf leather used for the inner lining ensure the highest wearing comfort.	
Shoe Designs	REGULAR left and right SLIM left and right	REGULAR left/SLIM right SLIM left/REGULAR right



## PARKSTREET

with Heelless Outsole Made of EVA Material  
with Rubber Coating

## URBANSTREET

with Heelless Outsole Made of EVA Material  
with Rubber Coating

Shoe Measurements		Use
pitch*	0mm	h-Cast
heel height**	14mm	Orthosis Configurator

Size	32-48	
Material	Sturdy fabric made of densely woven cotton – durable, easy to clean, guarantees a carefree use in everyday life.	



## CITYDRIVE

with Light, Durable Outsole Made of Polyurethan (PU)

Shoe Measurements		Use
pitch*	10mm	h-Cast
heel height**	25mm	Orthosis Configurator

Size	36-48	
Material	Fine Nubuck leather used for the upper material and soft calf leather for the inner lining ensure the highest wearing comfort.	

\* The height of the pitch corresponds to the height of the h-Cast plate to be used in the physical examination.

\*\* The heel height is used in the Orthosis Configurator as part of the calculation basis.

All orthosis shoes are sold in pairs for each shoe size.