



CROSSROADS

Orthosis shoe. Different.





CROSSROADS

Style and function for everyday life.

The CROSSROADS orthosis shoe is different. It is the result of years of experience, which we gained during the construction of custom-made orthoses. The requirements and demands for an orthosis shoe of people wearing orthoses were the main focus.



A unisex shoe that stands out from common orthosis shoes not only because of its extraordinary design features but as well because of its modern style. Casual through everyday life with CROSSROADS – the original with the boomerang.

CROSSROADS. Handmade with compelling details.



What Is Different?
The integrated elastic textile material offers comfortable space for orthoses with ankle joint and at the same time prevents the deformation of the shoe's upper.



Stability
The higher, anatomically formed upper and the integrated heel cap guarantee high stability.

Heel Strap
The heel strap facilitates putting on and taking off the shoe.

Guide for One-Handed Lacing
Each pair of shoes comes with a Guide for One-Handed Lacing, which offers support for patients with hemiplegia when lacing their shoes.

Cow Leather
Calf Leather

High Quality Leather
Certified cow leather used for the upper material and calf leather used for the inner lining ensure the highest wearing comfort.



2 Insoles
The two insoles (5mm each) are removable and serve as volume compensation. The upper sole is leather-coated for a pleasant foot climate.

Multifunctional Outsole
The outsole, which is made of EVA material, insulates from coldness and wetness and has an elastic anterior third as well as a slip-resistant, shock-absorbing tread.

CROSSROADS. Comfortable fit with and without orthosis.

CROSSROADS orthosis shoes always guarantee a perfect fit thanks to their special design and the carefully chosen material. The integrated elastic textile material as well as two removable insoles offer a flexible volume. Furthermore, it is possible to create an outsole elevation with an EVA sheet or to attach a vibram® flat sole under the shoe for some extra slip resistance. CROSSROADS is available as a pair in four different combinations of the widths REGULAR and SLIM, which makes the fit ideal for any foot – with and without orthosis.



Order Information



REGULAR left and right

for normal feet with
and without orthosis

Art. No.	Insole Length
OS4532-36	235mm
OS4532-37	242mm
OS4532-38	249mm
OS4532-39	254mm
OS4532-40	261mm
OS4532-41	268mm
OS4532-42	275mm
OS4532-43	281mm
OS4532-44	287mm
OS4532-45	294mm
OS4532-46	301mm
OS4532-47	307mm
OS4532-48	314mm



SLIM | REGULAR left | right

for slender feet with
orthosis on the right leg

Art. No.	Insole Length
OS8932-36	235mm
OS8932-37	242mm
OS8932-38	249mm
OS8932-39	254mm
OS8932-40	261mm
OS8932-41	268mm
OS8932-42	275mm
OS8932-43	281mm
OS8932-44	287mm
OS8932-45	294mm
OS8932-46	301mm
OS8932-47	307mm
OS8932-48	314mm



REGULAR | SLIM left | right

for slender feet with
orthosis on the left leg

Art. No.	Insole Length
OS9932-36	235mm
OS9932-37	242mm
OS9932-38	249mm
OS9932-39	254mm
OS9932-40	261mm
OS9932-41	268mm
OS9932-42	275mm
OS9932-43	281mm
OS9932-44	287mm
OS9932-45	294mm
OS9932-46	301mm
OS9932-47	307mm
OS9932-48	314mm



SLIM left and right

for slender feet with
orthosis without joints

Art. No.	Insole Length
OS4132-36	235mm
OS4132-37	242mm
OS4132-38	249mm
OS4132-39	254mm
OS4132-40	261mm
OS4132-41	268mm
OS4132-42	275mm
OS4132-43	281mm
OS4132-44	287mm
OS4132-45	294mm
OS4132-46	301mm
OS4132-47	307mm
OS4132-48	314mm

vibram® Flat Sole



Shoe Size	Art. No.
37/38	FS9701-3738
39/40	FS9701-3940
41/42	FS9701-4142
43/44	FS9701-4344
45/46	FS9701-4546
47/48	FS9701-4748

EVA Sheet for Outsole Elevation

Art. No.
PL4532-0S
Length x Width x Thickness
400 x 200 x 20mm

Pitch and Heel Height

Shoe Measurements	Use
pitch*	15mm h-Cast
heel height**	25mm Orthosis Configurator

* The height of the pitch corresponds to the height of the h-Cast plate to be used in the physical examination.

** The heel height is used in the Orthosis Configurator as part of the calculation basis.