


NEURO CLASSIC-SPRING

System Ankle Joint with Conversion Option to plug + go Modularity



available in steel and titanium

can be converted into a system joint
 with plug + go modularity
 by exchanging the functional unit
 and filing the system stirrup

improved integration into
 the laminate thanks to the
 slanted spring duct

exceptionally long lifespan
 due to redesign
 of the pressure springs



straight



bent inwards proximally



bent outwards proximally



10 mm



12 mm



14 mm



16 mm



20 mm


NEURO VARIO-SPRING 2

System Ankle Joint with plug + go Modularity

System joints with **plug + go modularity** are equipped with identical system stirrups, joint's upper parts and assembly/lamination dummies. All functional differences derive from the cover plate, which can be exchanged quickly and easily to convert between the system joints. System joints with **plug + go modularity** include the **NEURO VARIO 2**, the **NEURO VARIO-SWING** and the **NEURO SWING 2**.



NEURO VARIO-SPRING 2



improved integration into the laminate thanks to the slanted spring duct

integrated noise reduction

exceptionally long lifespan due to redesign of the pressure springs

available in steel and titanium



straight



bent inwards proximally



bent outwards proximally



10 mm



12 mm



14 mm



16 mm



20 mm