

Adjusting spring force of ankle foot orthoses according to gait type helps improving joint kinematics and time-distance parameters in patients with hemiplegia following stroke

Daniel Sabbagh, Jörg Fior, Ralf Gentz

Adjusting spring force of ankle foot orthoses according to gait type helps improving joint kinematics and time-distance parameters in patients with hemiplegia following stroke



Stroke

impairment of locomotion programmes

malfunction of executing extremities

false biomechanical situation

compensatory mechanisms

pathological gait

accompanying spasticity


ankle-foot orthoses (AFOs)

Gait classification in stroke

CLASSIFICATION OF GAIT TYPES ACCORDING TO PERRY					
Class	Description	Velocity		Knee (MSt)	Ankle (MSw)
Class I	FAST Walker	Fast	(55%)	Normal	Neutral position
Class II	MODERATE Walker	Moderate	(32%)	Flexed	Neutral position
Class III	SLOW EXTENDED Walker	Slow	(17%)	Hyperextended	Plantar flexion
Class IV	SLOW FLEXED Walker	Slow	(9%)	Flexed	Plantar flexion

[Perry et al. 2014]

N.A.P.[®] Gait Classification



[Sabbagh et al. 2014]

Orthotic management in stroke rehabilitation

AFOs play an important role in stroke management [Fatone 2009]

different AFO designs [Sabbagh et al. 2013]

different outcome parameters

several studies on orthotic intervention in stroke [Bowers et al. 2004, Condie et al. 2008]

- positive effects of AFOs – but not throughout
- mostly no differentiation of gait types
- poor methodology (e. g. different AFO designs)

positive results of changing AFOs spring force [Kobayashi et al. 2012, Kerkum et al. 2015]

Central Question: Can gait be improved by setting the AFOs stiffness according to the gait type?

Patients

n=8 (ø age 52.4, ø weight 82 kg, ø height 177 cm)

ischemic insult (Middle Cerebral Artery Stroke)

hemiplegia

gait type 1a+b (n=5), gait type 2a+b (n=3)

Inclusion criteria

6 Minute Walking Test, TUG Test

Exclusion criteria

pain

walking aids

Gait analysis

2-dimensional video analysis

2 conditions: a) Standardised footwear + DA-AFO, b) shoes only



3 full gait cycles

- time-distance parameters
- lateral kinematics: hip, knee, ankle
- maximum joint positions in stance (0-65% of gait cycle)

Wilcoxon rank-sum test (* $p < 0.05$, ** $p < 0.01$)

DA-AFO

Gait type 1 (knee hyperextension)




Gait type 2 (knee hyperflexion)

Kinematics

Gait type 1 (n=5)

Gait type 2 (n=3)



Stance: Flexion +

2nd half of stance: Flexion -

Stance: Flexion +

2nd half of stance: Flexion -

1st rocker improved
Stance: DF +
Swing: DF +


1st rocker improved
2nd half of stance: DF -
Terminal stance: Peak DF -

Time-Distance Parameters

			Gait type 1 (n=5)			Gait type 2 (n=3)			Reference (n=26)
			Shoes only	DA-AFO	W RST	Shoes only	DA-AFO	W RST	Shoes only
Hip	max. ext	[°]	-10.2 (±6.1)	-9.9 (±3.6)	-	-13.1 (±0.9)	-15.2 (±1.7)	-	-19.0 (±4.5)
Knee	max. ext	[°]	0.8 (±4.0)	5.2 (±3.2)	*	11.4 (±3.3)	6.4 (±1.2)	*	3.8 (±4.1)
Ankle	at IC	[°]	-9.9 (±5.0)	4.6 (±2.1)	**	-5.1 (±11.4)	1.2 (±4.6)	-	-1.2 (±3.8)
	max. DF	[°]	10.7 (±3.1)	11.7 (±2.8)	-	18.1 (±1.9)	13.9 (±1.3)	*	14.7 (±4.3)
Stride length		[m]	0.9 (±0.3)	1.0 (±0.2)	*	0.8 (±0.3)	0.9 (±0.3)	*	1.6 (±0.1)
Velocity		[m/s]	0.5 (±0.3)	0.7 (±0.2)	**	0.5 (±0.2)	0.5 (±0.2)	*	1.4 (±0.2)
Cadence		[Steps/min]	65.6 ±15.1	73.7 (±8.7)	*	72.6 (±4.3)	74.7 (±8.7)	-	104.4 (±9.0)
Stance		[% GC]	67.5 (±8.5)	66.0 (±7.1)	-	70.6 (±5.4)	64.6 (±4.0)	*	65.3 (±2.0)
Swing		[% GC]	32.5 (±8.5)	34.0 (±7.1)	-	29.4 (±5.4)	35.5 (±4.0)	*	34.7 (±2.0)

Wilcoxon rank-sum test (W RST): * indicates $p < 0.05$, ** indicates $p < 0.01$


Heel Rocker




Dorsal spring

Gait type 1

Gait type 2




Dorsiflexion resistance



Ventral spring

Gait type 1

Gait type 2



Conclusions

improvements in both gait types

time-distance parameters follow kinematic improvements

different spring forces different effects individual biomechanical situations

setting spring force according to gait type leads to gait improvements [Kerkum et al. 2015]

individual setting and alignment of AFO is necessary



**Thank you
for your
attention!**



Patient #	Gait type	Spring unit					
		ventral			dorsal		
		Code	Spring force	Nm/deg	Code	Spring force	Nm/deg
1	1	yellow	very strong	0.64	yellow	very strong	0.64
2	1	yellow	very strong	0.64	yellow	very strong	0.64
3	2	yellow	very strong	0.64	green	medium	0.24
4	2	white	strong	0.52	green	medium	0.24
5	1	red	extra strong	1.71	green	medium	0.24
6	1	white	strong	0.52	green	medium	0.24
7	2	red	extra strong	1.71	green	medium	0.24
8	1	yellow	very strong	0.64	yellow	very strong	0.64